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Peter Coleman  
Sr Reg Segment Marketing Manager  
Asia Pacific Japan, Illumina

© 2016 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiniSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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


1. Minasi G M, *et al*; Correlation between aneuploidy, standard morphology evaluation and morphokinetic development in 1730 biopsied blastocysts: a consecutive case series study. *Human Reproduction*, Volume 31, Issue 10, 1 October 2016, Pages 2245–2254,
2. Machtinger R, Racowsky C. Morphological systems of human embryo assessment and clinical evidence. *Reproductive BioMedicine Online* (2013) 26, 210– 221
3. Kovacs P. Embryo selection: the role of time-lapse monitoring. *Reproductive Biology and Endocrinology* 2014, 12:124

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2008

2014

# Diminished effect of maternal age on implantation after preimplantation genetic diagnosis with array comparative genomic hybridization

Gary L. Harton, B.S.,<sup>a</sup> Santiago Munné, Ph.D.,<sup>b</sup> Mark Surrey, M.D.,<sup>c</sup> Jamie Grifo, M.D., Ph.D.,<sup>d</sup> Brian Kaplan, M.D.,<sup>e</sup> David H. McCulloh, Ph.D., H.C.L.D.,<sup>d</sup> Darren K. Griffin, Ph.D.,<sup>f</sup> and Dagan Wells, Ph.D.,<sup>g,h</sup> for the PGD Practitioners Group<sup>i</sup>

<sup>a</sup> Bluegnome, La Jolla, California; <sup>b</sup> Reprogenetics, Livingston, New Jersey; <sup>c</sup> Southern California Reproductive Center, Beverly Hills, California; <sup>d</sup> NYU Fertility Center, New York, New York; <sup>e</sup> Highland Park IVF Center, Fertility Centers of Illinois, Highland Park, Illinois; <sup>f</sup> School of Biosciences, University of Kent, Canterbury, United Kingdom; and <sup>g</sup> Reprogenetics UK and <sup>h</sup> Nuffield Department of Obstetrics and Gynaecology, University of Oxford, Oxford, United Kingdom; and <sup>i</sup> Centers in the PGD Practitioners Group are listed at the end of the article

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
1. Harton GL, Munne S, Surrey M, et al. Diminished effect of maternal age on implantation after preimplantation genetic diagnosis with array comparative genomic hybridization. *Fertil Steril*. 2013;100(6):1695-1703.

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\*\*Statistical significance of difference in rates between age groups - Aneuploidy  $p < 0.001$ ; Implantation  $p > 0.05$ ; Miscarriage  $p > 0.25$

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1. Harton GL, Munne S, Surrey M, et al. Diminished effect of maternal age on implantation after preimplantation genetic diagnosis with array comparative genomic hybridization. *Fertil Steril.* 2013;100(6):1695-1703.

## Single thawed euploid embryo transfer improves IVF pregnancy, miscarriage, and multiple gestation outcomes and has similar implantation rates as egg donation

Jamie A. Grifo • Brooke Hodes-Wertz • Hsiao-Ling Lee •  
Esmeralda Amperloquio • Melicia Clarke-Williams • Alexis Adler

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
- Statistically significant

1. Grifo JA, Hodes-Wertz B, Lee HL, Amperloquio E, Clarke-Williams M, Adler A. Single thawed euploid embryo transfer improves IVF pregnancy, miscarriage, and multiple gestation outcomes and has similar implantation rates as egg donation. *J. Assist. Reprod. Genet.* 2013; 30, 259–64.

## REPRODUCTIVE ENDOCRINOLOGY AND INFERTILITY

# Obstetrical and neonatal outcomes from the BEST Trial: single embryo transfer with aneuploidy screening improves outcomes after in vitro fertilization without compromising delivery rates

Eric J. Forman, MD; Kathleen H. Hong, MD; Jason M. Franasiak, MD; Richard T. Scott Jr, MD

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
1. Forman EJ, Hong KH, Franasiak JM, Scott RT. Obstetrical and neonatal outcomes from the BEST Trial: single embryo transfer with aneuploidy screening improves outcomes after in vitro fertilization without compromising delivery rates. *Am. J. Obstet. Gynecol.* 2014; 210, 157.e1–6.

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## Development and validation of a next-generation sequencing–based protocol for 24-chromosome aneuploidy screening of embryos

Francesco Fiorentino, Ph.D.,<sup>a</sup> Anil Biridik, M.Sc.,<sup>a</sup> Sara Bono, B.Sc.,<sup>a</sup> Letizia Spizzichino, B.Sc.,<sup>a</sup> Ettore Cotroneo, B.Sc.,<sup>a</sup> Giuliano Cottone, B.Sc.,<sup>a</sup> Felix Kokocinski, Ph.D.,<sup>b</sup> and Claude-Edouard Michel, Ph.D.<sup>b</sup>

<sup>a</sup> Genoma Molecular Genetics Laboratory, Rome, Italy; and <sup>b</sup> Bluegenome, Cambridge, United Kingdom

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1. Fiorentino F, Biridik A, Bono S, et al Development and validation of a next-generation sequencing-based protocol for 24-chromosome aneuploidy screening of embryos. 2014; 101(5):1375-82.

## Application of next-generation sequencing technology for comprehensive aneuploidy screening of blastocysts in clinical preimplantation genetic screening cycles

Francesco Fiorentino<sup>1,\*</sup>, Sara Bono<sup>1</sup>, Anil Biricik<sup>1</sup>, Andrea Nuccitelli<sup>1</sup>, Ettore Cotroneo<sup>1</sup>, Giuliano Cottone<sup>1</sup>, Felix Kokocinski<sup>2</sup>, Claude-Edouard Michel<sup>2</sup>, Maria Giulia Minasi<sup>3</sup>, and Ermanno Greco<sup>3</sup>

<sup>1</sup>Molecular Genetics Laboratory, 'GENOMA', Via di Castel Giubileo, 11, 00138 Rome, Italy <sup>2</sup>Illumina, Inc., Cambridge, UK  
<sup>3</sup>Reproductive Medicine, European Hospital, Via Portuense, 700, 00149 Rome, Italy

\*Correspondence address. Molecular Genetics Laboratory, 'GENOMA', Via di Castel Giubileo, 11, 00138 Rome, Italy.  
E-mail: fiorentino@laboratoriotgenoma.it

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1. Fiorentino F, Bono S, Biricik A, et al. Application of next-generation sequencing technology for comprehensive aneuploidy screening of blastocysts in clinical preimplantation genetic screening. 2014; 29(12):2802-13.

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

\* Per embryo transferred

1. Harton G, Munné S, Surrey M, et al. Diminished effect of maternal age on implantation after Preimplantation Genetic Diagnosis with array comparative genomic hybridization. *Fertil Steril.* 2013; 100(6):1695-703.
2. Grifo JA, Hodes-Wertz B, Lee HL, Amperloquio E, Clarke-Williams M, Adler A. Single thawed euploid embryo transfer improves IVF pregnancy, miscarriage, and multiple gestation outcomes and has similar implantation rates as egg donation. *J. Assist. Reprod. Genet.* 2013; 30, 259–64.
3. Scott RT, Upham KM, Forman EF, et al. Blastocyst biopsy with comprehensive chromosome screening and fresh embryo transfer significantly increases IVF IRs and delivery rates: a randomized controlled trial. *Fertil. Steril.* 2013; 100(3):697-703
4. Forman EJ, Hong KH, Ferry KM, et al. In vitro fertilization with single euploid blastocyst transfer: a randomized controlled trial. *Fertil Steril.* 2013; 100(1):100-7.
5. Yang Z, Liu J, Collins GC, et al. Selection of single blastocysts for fresh transfer via standard morphology assessment alone and with array CGH for good prognosis IVF patients: results from a randomized pilot study. *Mol. Cytogenet.* 2012; 5(1):24

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American College of Obstetricians and Gynecologists. Screening for fetal aneuploidy. Practice Bulletin No. 163. *Obstet Gynecol.* 2016; 127(5):e123-37.

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
1. Barrett, A, Zimmerman BG, Wang D, Holloway A, Chitty L. Implementing prenatal diagnosis based on cell-free fetal DNA: Accurate identification of factors affecting fetal DNA yield. *PLoS One*. 2011;6(10):e25202..
2. Nigam A, Saxena P, Prakash A, Acharya A. Detection of fetal nucleic acid in maternal plasma: A novel noninvasive prenatal diagnostic technique. *J Internl Med Sci Acad*. 2012; 25(3): 119-120.

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

1. Taneja PA, Snyder HL, de Feo E, et al. Noninvasive prenatal testing in the general obstetric population: clinical performance and counseling considerations in over 85 000 cases. *Prenat Diagn.* 2016; 36(3):237-243.

2. McCullough RM, Almasri EA, Guan X, et al. Non-invasive prenatal chromosomal aneuploidy testing – clinical experience: 100 000 clinical samples. *PLoS One.* 2014; 9(10):e109173.

3. Juneau K, Mohseni M, Wang ET, et al. Microarray-based cell-free DNA analysis improves noninvasive prenatal testing. *Fetal Diagn Ther.* 2014; 36(4):282-286.

4. Stokowski R, Wang E, White K, et al. Clinical performance of non-invasive prenatal testing (NIPT) using targeted cell free DNA analysis in maternal plasma with microarrays or next generation sequencing (NGS) is consistent across multiple controlled clinical studies. *Prenat Diagn.* 2015; 35(12):1243-1246.

5. White K, Wang E, Batey A, et al. Performance of targeted cell-free DNA analysis with microarray quantitation for assessment of fetal sex and sex chromosome aneuploidy risk.

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
1. Norton ME, Jacobsson B, Swamy GK, et al. Cell-free DNA analysis for noninvasive examination of trisomy. *New Engl J Med*. 2015; 372(17):1589-1597.

2. Ryan A, Hunkapiller N, Banjevic M, et al. Validation of an enhanced version of a single -nucleotide polymorphism-based noninvasive prenatal test for detection of fetal aneuploidies. *Fetal Diagn Ther*. 2016;doi:10.1159/000442931.

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
- ACMG 2016

American College of Obstetricians and Gynecologists. Screening for fetal aneuploidy. Practice Bulletin No. 163. *Obstet Gynecol.* 2016; 127(5):e123-137.

Benn P, Borrell A, Chiu RWK, et al. Position statement from the Chromosome Abnormality Screening Committee on behalf of the Board of the International Society for Prenatal Diagnosis. *Prenat Diagn.* 2015;35(8):725-734. doi:10.1002/pd.4608.


Gregg AR, Skotko BG, Benkendorf JL, et al. Noninvasive prenatal screening for fetal aneuploidy, 2016 update: a position statement of the American College of Medical Genetics and Genomics. *Genet Med.* 2016: doi:10.1038/gim.2016.97.

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
Data compiled from published clinical studies. WGS Illumina<sup>1</sup>, WGS other<sup>2</sup>, Targeted SNP Sequencing.<sup>3</sup>

\*There are no published studies assessing failure rate for NIPT assays performed using SNP arrays in clinical samples.

1. Taneja P, Snyder H, de Feo E, et al. Noninvasive prenatal testing in the general obstetric population: clinical performance and counseling considerations in over 85,000 cases. *Prenat Diagn.* 2016;36:237-43.
2. McCullough RM, Almasri EA, Guan X, et al. Non-invasive prenatal chromosomal aneuploidy testing--clinical experience: 100,000 clinical samples. *PLoS One.* 2014;9(10):e109173.
3. Ryan A, Hunkapiller N, Banjevic M, et al. Validation of an enhanced version of a single-nucleotide polymorphism-based noninvasive prenatal test for detection of fetal aneuploidies. *Fetal Diagn Ther.* 2016;doi:10.1159/000442931.

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Lên tới


**65%**

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1. Pergament E, Cuckle H, Zimmermann B, et al. Single-nucleotide polymorphism-based noninvasive prenatal screening in a high-risk and low-risk cohort. *Obstet Gynecol.* 2014; 124(2 Pt 1):210-218.

2. American College of Obstetricians and Gynecologists. Screening for fetal aneuploidy. Practice Bulletin No. 163. *Obstet Gynecol.* 2016; 127(5):e123-37.

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Test Failures May Lead to Invasive Procedures.†—Theoretical example of the number of invasive procedures requested due to NIPT failure. Failure rates include assay failures and samples rejected due to low fetal fraction. Assay failure rate for the Harmony test is based on NGS studies and may not be consistent with actual test results achieved using the array-based Harmony Test currently in use.

† Affected pregnancies with a screening test failure were excluded from the number of detected T21.

## References

3. Taneja P, Snyder H, de Feo E, et al. Noninvasive prenatal testing in the general obstetric population: clinical performance and counseling considerations in over 85,000 cases. *Prenat Diagn*. 2016;36:237-43.
4. Ryan A, Hunkapiller N, Banjevic M, et al. Validation of an enhanced version of a single-nucleotide polymorphism-based noninvasive prenatal test for detection of fetal aneuploidies. *Fetal Diagn Ther*. 2016;doi:10.1159/000442931.
5. McCullough RM, Almasri EA, Guan X, et al. Non-invasive prenatal chromosomal aneuploidy testing—clinical experience: 100,000 clinical samples. *PLoS One*. 2014;9(10):e109173.

12. Norton ME, Jacobsson B, Swamy GK, et al. Cell-free DNA analysis for noninvasive examination of trisomy. *N Engl J Med*. 2015;doi: 10.1056/NEJMoa1407349.

18. Nicolaidis KH. Screening for fetal aneuploidies at 11 to 13 weeks. *Prenat Diagn*. 2011;31(1):7-15.

19. Wald NJ, Rodeck C, Hackshaw AK, et al. SURUSS in perspective. *Semin Perinatol*. 2005;29(4):225-235.

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1. Taneja PA, Snyder HL, de Feo E, et al. Noninvasive prenatal testing in the general obstetric population: clinical performance and counseling considerations in over 85 000 cases. *Prenat Diagn.* 2016; 36(3):237-243.
2. Rava RP, Srinivasan A, Sehnert AJ, Bianchi DW. Circulating fetal cell-free DNA fractions differ in autosomal aneuploidies and monosomy X. *Clin Chem.* 2014; 60(1):243-50.
3. Gil MM, Quezada MS, Revello R, Akolekar R, Nicolaides KH. Analysis of cell-free DNA in maternal blood in screening for fetal aneuploidies: updated meta-analysis. *Ultrasound Obstet Gynecol.* 2015; 45(3):249-266.

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Câu hỏi?

[pcoleman@illumina.com](mailto:pcoleman@illumina.com)